

For Research Use Only

IL-12A/IL-12 p35 Recombinant antibody, PBS Only (Capture)

Catalog Number: 85754-3-PBS



Basic Information

Catalog Number:

85754-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC104982

GeneID (NCBI):

3592

UNIPROT ID:

P29459

Full Name:

interleukin 12A (natural killer cell stimulatory factor 1, cytotoxic lymphocyte maturation factor 1, p35)

Calculated MW:

25 kDa

Purification Method:

Protein A purification

CloneNo.:

242887B7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Cytometric bead array standard curve of MP02056-1, IL-12A/IL-12 p35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85754-3-PBS. Detection antibody: 85754-2-PBS. Standard: Eg3341. Range: 0.781-50 ng/mL



Cytometric bead array standard curve of MP02056-2, IL-12A/IL-12 p35 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 85754-3-PBS. Detection antibody: 85754-1-PBS. Standard: Eg3341. Range: 0.781-100 ng/mL