

For Research Use Only

CD247 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86659-2-PBS



Basic Information

Catalog Number: 86659-2-PBS	GenBank Accession Number: BC025703	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 919	CloneNo.: 251417H4
Source: Rabbit	UNIPROT ID: P20963	
Isotype: IgG	Full Name: CD247 molecule	
Immunogen Catalog Number: EG5018	Calculated MW: 16 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD247 concentration was determined to be 99.2 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02576-1, CD247 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86659-2-PBS. Detection antibody: 86659-1-PBS. Standard: Eg5018. Range: 0.469-60 ng/mL.