

For Research Use Only

# Anti-Human TNFR2/CD120b Rabbit Recombinant Antibody



Catalog Number: 98010-2-RR

## Basic Information

|                                      |   |   |
|--------------------------------------|---|---|
| <b>Catalog Number:</b><br>98010-2-RR | <b>GenBank Accession Number:</b><br>BC052977                                  | <b>Purification Method:</b><br>Protein A purification |
| <b>Size:</b><br>100ug, 1000 µg/ml    | <b>GeneID (NCBI):</b><br>7133   | <b>CloneNo.:</b><br>230328B6                          |
| <b>Source:</b><br>Rabbit             | <b>UNIPROT ID:</b><br>P20333  |   |
| <b>Isotype:</b><br>IgG               | <b>Full Name:</b><br>tumor necrosis factor receptor<br>superfamily, member 1B |   |
|                                      | <b>Calculated MW:</b><br>48 kDa   |   |

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Human

## Background Information

Tumor necrosis factor-alpha (TNFA/TNFSF2) is a multifunctional cytokine that plays a key role in regulating inflammation, immune functions, host defense, and apoptosis (PMID: 16407280). TNFA signals through two distinct cell surface receptors, TNFR1 (TNFRSF1A, CD120a, p55) and TNFR2 (TNFRSF1B, CD120b, p75). TNFR1 is widely expressed, whereas TNFR2 exhibits more restricted expression, being found on CD4 and CD8 T lymphocytes, endothelial cells, microglia, oligodendrocytes, neuron subtypes, cardiac myocytes, thymocytes and human mesenchymal stem cells (PMID: 20489699; 22374304). In contrast to TNFR1, TNFR2 does not have a death domain. TNFR2 only signals for antiapoptotic reactions. However, recent evidence indicates that TNFR2 also signals to induce TRAF2 degradation (PMID: 22374304). Various defects in the TNFR2 pathway, due to polymorphisms in the TNFR2 gene, upregulated expression of TNFR2 and TNFR2 shedding, have been implicated in the pathology of several autoimmune disorders (PMID: 20489699).

## Storage

**Storage:**  
Store at 2-8°C. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



$1 \times 10^6$  human peripheral blood leukocytes were surface stained with 0.25  $\mu$ g Anti-Human TNFR2/CD120b Rabbit Recombinant Antibody (98010-2-RR, Clone:230328B6) and PE-conjugated Goat Anti-Rabbit IgG, or 0.25  $\mu$ g Isotype Control. Cells were not fixed.



Surface Plasmon Resonance (SPR) kinetic assays of 98010-2-RR against Human TNFR2/CD120b were performed. The affinity constant is 0.4046 nM.