

For Research Use Only

Anti-Human uPAR/CD87 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98048-1-PBS



Basic Information

Catalog Number:

98048-1-PBS

Size:

1mg, 2 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_002659.4

GeneID (NCBI):

5329

UNIPROT ID:

Q03405-1

Full Name:

PLAUR

Calculated MW:

37 kDa

Purification Method:

Protein A purification

CloneNo.:

240366G10

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 untreated (black) or PMA-treated U-937 cells were surface stained with 0.25 μg Anti-Human uPAR/CD87 Rabbit Recombinant Antibody (98048-1-RR, Clone: 240366G10) and PE-conjugated Goat Anti-Rabbit IgG (red). 1×10^6 PMA treated U-937 cells were surface stained with 0.25 μg Isotype Control (blue). Cells were not fixed. This data was developed using the same antibody clone with 98048-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98048-1-RR against Human uPAR were performed. The affinity constant is 4.70 nM.