

For Research Use Only

# Anti-Human VE-cadherin/CD144 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98071-1-PBS



## Basic Information

<b>Catalog Number:</b> 98071-1-PBS	<b>GenBank Accession Number:</b> NM_001795.5	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1mg, 2mg/ml	<b>GeneID (NCBI):</b> 1003	<b>CloneNo.:</b> 240755B2
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P33151-1	
<b>Isotype:</b> IgG	<b>Full Name:</b> cadherin 5, type 2 (vascular endothelium)	
	<b>Calculated MW:</b> 88 kDa	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
human

## Background Information

Cadherins are a family of transmembrane glycoproteins that mediate calcium-dependent cell-cell adhesion and play an important role in the maintenance of normal tissue architecture. Vascular endothelial cadherin (VE-cadherin), also known as Cadherin-5 (CDH5) or CD144, is a member of the type II classical cadherin family of cell adhesion proteins (PMID: 21269602). VE-cadherin is expressed specifically in endothelial cells and mediates homophilic adhesion in the vascular endothelium (PMID: 1522121; 8555485; 21269602). VE-cadherin plays a role in the organization of lateral endothelial junctions and in the control of permeability properties of vascular endothelium (PMID: 1522121). VE-cadherin has also been shown to be required for angiogenesis (PMID: 16473763; 18162609). The calculated molecular weight of VE-cadherin is 88 kDa and the apparent molecular weight of 120-140 kDa is higher due to post-translational glycosylation and phosphorylation (PMID: 10460833; 29894844). Full-length VE-cadherin can be proteolytically cleaved to generate a fragment of 90-100 kDa (PMID: 9786462; 22064597).

## Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



$1 \times 10^6$  HUVEC cells were stained with 0.25  $\mu$ g Anti-Human VE-cadherin/CD144 Rabbit Recombinant Antibody (98071-1-RR, Clone:240755B2) (red), or 0.25  $\mu$ g Isotype Control (blue). Cells were not fixed. This data was developed using the same antibody clone with 98071-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98071-1-RR against Human VE-cadherin/CD144 were performed. The affinity constant is below 1 pM.