

For Research Use Only

Anti-Mouse CD27 Rabbit Recombinant Antibody, PBS Only

Catalog Number: 98160-1-PBS



Basic Information

Catalog Number: 98160-1-PBS	GenBank Accession Number: NM_001033126.2	Purification Method: Protein A purification
Size: 1mg, 2 mg/ml	GeneID (NCBI): 21940	CloneNo.: 241378E11
Source: Rabbit	UNIPROT ID: P41272	
Isotype: IgG	Full Name: CD27 antigen	
	Calculated MW: 28kDa	

Applications

Tested Applications:
FC

Species Specificity:
mouse

Background Information

CD27 (also known as TNFRSF7) is a type I glycoprotein expressed on some B cells and the majority of T cells, and is a member of the tumor necrosis factor (TNF) receptor family. CD27 is required for generation and long-term maintenance of T cell immunity (PMID: 11062504). It is a receptor for CD70 (CD27L). Ligation of CD27 by CD70 induces strong ubiquitination of TRAF and the activation of both canonical and non-canonical NF-kappaB pathways, as well as the JNK pathway (PMID: 19426224). CD27 may also play a role in apoptosis through association with SIVA1.

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS Only

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD27 Rabbit Recombinant Antibody (98160-1-RR, Clone:241378E11) or Isotype Control, and PE-Conjugated Goat Anti-Rabbit IgG(H+L), cells were co-stained with APC Anti-mouse CD3. Cells were not fixed. This data was developed using the same antibody clone with 98160-1-PBS in a different storage buffer formulation.



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Anti-Mouse CD27 Rabbit Recombinant Antibody (98160-1-RR, Clone: 241378E11) (red), or Isotype Control (blue), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 98160-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 98160-1-RR against Mouse CD27 were performed. The affinity constant is 1.87 nM.