

For Research Use Only

Human IgG1, κ isotype control Recombinant Antibody, PBS Only

Catalog Number: 99026-1-PBS



Basic Information

Catalog Number:

99026-1-PBS

Concentration:

1mg, 3 mg/ml

Source:

Human

Isotype:

IgG1

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

262113A1

Applications

Tested Applications:

FC

Species Specificity:

human

Background Information

This product is an isotype control antibody, which is used to estimate the non-specific binding of target primary antibodies due to Fc receptor binding or other protein-protein interactions. An isotype control antibody should have the same immunoglobulin type and be used at the same concentration as the test antibody.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ mouse splenocytes were surface stained with 0.25 ug Human IgG1 kappa isotype control RecAb (99026-1-HR, Clone: 262113A1) (red in histogram), 0.25 ug Anti-Mouse CD23 Rabbit Recombinant Antibody (98164-1-RR, Clone: 241385D1) (green in histogram), or unstained cells (blue in histogram), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were not fixed. This data was developed using the same antibody clone with 99026-1-PBS in a



1x10⁶ human PBMCs were surface stained with 0.25 ug Human IgG1 kappa isotype control RecAb (99026-1-HR, Clone: 262113A1) (red in histogram), or 0.25 ug Anti-Human CD4 Rabbit Recombinant Antibody (98042-1-RR, Clone: 240427E12) (green in histogram), or unstained (blue in histogram), and PE-Conjugated Goat Anti-Rabbit IgG(H+L). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed. This data was developed

