

For Research Use Only

# APC Anti-Human CD42b (AK2)

Catalog Number: APC-65163



## Basic Information

**Catalog Number:**

APC-65163

**Size:**

100tests, 5 µl/test

**Source:**

Mouse

**Isotype:**

IgG1

**GenBank Accession Number:**

BC027955

**GeneID (NCBI):**

2811

**ENSEMBL Gene ID:**

ENSG00000185245

**UNIPROT ID:**

P07359

**Full Name:**

glycoprotein Ib (platelet), alpha polypeptide

**Calculated MW:**

626 aa, 69 kDa

**Purification Method:**

Purified by protein-A affinity chromatography

**CloneNo.:**

AK2

**Excitation/Emission maximum wavelengths:**

650 nm / 660 nm

## Applications

**Tested Applications:**

FC

**Species Specificity:**

Human

## Background Information

CD42b, also known as platelet glycoprotein Ib alpha chain (GPIb alpha), is a type I transmembrane glycoprotein of 135-145 kDa (PMID: 2656709). It is expressed on platelets and megakaryocytes. CD42b and CD42c (GPIb beta) are linked by a disulphide bond to form GPIb (PMID: 7660135). GPIb forms a noncovalent complex with CD42a (GPIX) and CD42d (GPV) and acts as a receptor for von Willebrand factor (vWF) and thrombin (PMID: 7660135; 3759960). The GPIb-V-IX complex mediates vWF-dependent platelet adhesion to blood vessels. Defects in the expression of CD42b result in Bernard-Soulier syndromes and platelet-type von Willebrand disease (PMID: 9371310).

## Storage

**Storage:**

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> human peripheral blood platelets were surface co-stained with Anti-Human CD41a labeled with FlexAble CoraLite® Plus 488 Antibody Labeling Kit for Mouse IgG1 (KFA041) and 5 ul APC Anti-Human CD42b (APC-65163, Clone: AK2) or Mouse IgG1 Isotype Control. Cells were not fixed.