

For Research Use Only

APC Rabbit IgG Isotype Control Recombinant Antibody

Catalog Number: APC-98136



Basic Information

Catalog Number:

APC-98136

Size:

100ug, 200 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

240953C9

Excitation/Emission maxima
wavelengths:

650 nm / 660 nm

Applications

Tested Applications:

FC

Species Specificity:

n/a

Background Information

This product is an isotype control antibody, which is used to estimate the non-specific binding of target primary antibodies due to Fc receptor binding or other protein-protein interactions. An isotype control antibody should have the same immunoglobulin type and be used at the same concentration as the test antibody.

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ human peripheral blood leukocytes were unstained (left, blue in histogram) or surface stained with 0.1 ug APC Rabbit IgG Isotype Control Recombinant Antibody (APC-98136, Clone: 240953C9) (mid, red in histogram), or 0.1 ug APC Anti-Human CD45 (HI30) (APC-65109, Clone: HI30) (right, green in histogram). Cells were incubated with FC Receptor Block prior to staining. Cells were not fixed.