

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

RECOMBINANT HUMAN MMP9



Basic Information

Catalog Number:
Ag0552

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
E. coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Validated Application:
Blocking peptide

Peptide Sequence:

MSLWGQLVLLVLLVGGCFAAPRQRQSTLVLFPGDLRNLTDRL
AEEYLYRYGYTRVAEMRGESKSLGPALLLLQKQLSLPETGELDS
ATLKAMRTPRCG
(1-100 aa encoded by BC006093)

Selected Validation Data

For technical support and original validation data for this product, please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under
Proteintech Group brand and is not available to
purchase from any other manufacturer.