

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human TCEB3B protein



Basic Information

Catalog Number:

Ag10595

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

SATTALGDKQRKANESKGTRESWDSAKLPPVQESQ
SERLQAAGADSAGPKTVPNHVFSELWDLSEAWMQA
NYDPLSDSDSMTSQAKPEALSSPKFREEAAFPGRRVN
AKMPVYSGSRPACQLQVPTLRQQCAQVLRNNPDAL
SDVGEVPYWVLEPVLEGWRPDQLYRRKKDNHALVRE
TDELRRNHCFQDFKEEKQENKTWREQYLRLDAPE
QRLRVMTTNIRSARGNNPNGREAKMICFKSVAKTPYD
TSRRQEKSAGDADPENGEIKPASKPAGSSHTPSSQSS
SGGGRDSSSSILRWLPEKRANPCLSSSNEHAAPAAKT
RKQAAKKVAPLMAKAIRDYKRRFSRR
(402-753 aa encoded by BC036022)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.