

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# RECOMBINANT HUMAN TF

## Basic Information

**Catalog Number:**

Ag11668

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*E. coli*-derived, PET28a, with N-terminal 6\*His.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Validated Application:**

Blocking peptide

**Peptide Sequence:**

EGTCCQEAPTDECKPMKWCALSHHERLKCDEWSVNSVGKIECV  
SAETTEDCIAKIMNGEADAMSLDGGFVYIAGKCGLVPVLAENYN  
KSDNCEDTPEAGYFAVAWVKKSASDLTWDNLKGGKKSCHTAVGR  
TAGWNIPIMGLLYNKINHCRFDEFFSEGCAPGSKKDSCLKLKM  
GSGNLNLCFNNKEGYGYTGAFRCLVEKGDVAFVKHQTVFQN  
TGKKNPDPWAKNLNEKDYELLCLDGTTRKPMEEYANCHLARAP  
NHAVTRKDKKACVHKILRQQQLHFGSNVTDCCSGNFCLFRSET  
KDLLFRDDTVCLAKLHNRNTYKYLGEELYKAVGNLRKCSTSS  
LLEACTFRRP  
(347-698 aa encoded by BC059367)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under  
Proteintech Group brand and is not available to  
purchase from any other manufacturer.