

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# TLR3 Fusion Protein



## Basic Information

Catalog Number:

Ag12171

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Validated Application:

Blocking peptide

Peptide Sequence:

MRQTLPCIYFWGGLLPFGMLCASSTTKCTVSHEVAD  
CSHLKLTQVPDDLPTNITVLNLTHNQLRRLPAANFTRY  
SQLTSLDVGFNITISKLEPELCQKLPMLKVLNLQHNELS  
QLSDKTFAFCTNLTELHMSNSIQIKNNPFVKQKNLI  
TLDLSHNGLSSTKLGTVQVQLENLQELLSNNKIQALKS  
EELDIFANSSLKLELSSNQIKEFSPGCFHAI GRLFGLFL  
NNVQLGPSLTEKLCLELANTSIRNLSLSNSQLSTTSNTT  
FLGLKWTNLTMLDLSYNNLN VVGNDFAWLPQLEYFF  
LEYNNIQHLFSLHGLFNVRVYLNLRKRSFTKQSI SLASL  
PKIDD

(26-348 aa encoded by BC094737)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.