

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human DNAH6 protein



## Basic Information

Catalog Number:

Ag12774

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MTYTEQASKRHYCMRLTCFIRLNDYLIENTMHILTVNA  
VNSLLNHLTDKLRTPSADVIQKWITEEKPEVPDKKGT  
LMVEKQEEDESLLPMFLTELMLTVQSLLFEPSEDFLDG  
ILGAVNHCQNTVLSVPLVPDSYFDAFTSPYINNKLE  
GKTCGTGPSLAAVFEDDKNFHTIISQIKETIQAAFESA  
RIYAATFEKFQIFFKENESLDLQALKQEPDINFFSEQL  
EKYHKQHKDAVALRPTRNVGLLLIDTRLREKLIPLR  
CLEVLNFMPLRQSKKKVDAIIFEAQDAEYKLEFVPTTT  
TEYVHSLFLDEIQERIESELEDEGNIVTQMYKLEQYQ  
VPTPPEDFAVFATMKPSIVAVRNAIDKSVGDRESSIK  
QFCVHLGSDLEELNNEVNEVKLQAQVS  
(1-408 aa encoded by BC104884)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.