

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human Tyrosinase protein



## Basic Information

### Catalog Number:

Ag17314

### Size:

50 µg

### Form:

Available lyophilized

### Species:

human

### Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

### Biological Activity:

Not tested

### Endotoxin Level:

Please contact the lab for more information

### Peptide Sequence:

HFPRACVSSKNLMEKCCPPWSGDRSPCQLSGRGS  
CQNILLSNAPLGQFPFTGVDDRESWPSVFYNRTCQ  
CSGNFMGFNCGNCKFGFWGPNCTERRLLVRRNIFDL  
SAPEKDKFFAYLTAKHTISSDYVPIGTYGQMKNGST  
PMFNDINIYDLFVWMHYVSM DALLGGSEIWRDIDFA  
HEAPAFLPWHLRLLRWEQEIQKLTGDENFTIPYWDW  
RDAEKCDICTDEYMGQHPNPNLLSPASFFSSWQI  
VCSRLEEYNHQSCLNGTPEGPLRRNPGNHDKSRTPR  
LPSSADVEFCLSLTQYESGSM DKAANFSFRNTLE  
(19-345 aa encoded by BC027179)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.