

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human WFDC3 protein



Basic Information

Catalog Number:

Ag20207

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PET28a, with N-terminal 6*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GEHAKEGECPPHKNPCKELCQGDELCPAEQKCCTTG
CGRICRDI PKGRKRD C P R V I R K Q S C L K R C I T D E T C P G V
K K C C T L G C N K S C V V P I S K Q K L A E F G G E C P A D P L P C E E
L C D G D A S C P Q G H K C C S T G C G R T C L G D I E G G R G G D C P
K V L V G L C I V G C V M D E N C Q A G E K C C K S G C G R F C V P P
V L P P K L T M N P N W T V R S D S E I P V P
(25-231 aa encoded by BC148860)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.