

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human SLC9A3 protein



## Basic Information

**Catalog Number:**

Ag25637

**Size:**

50 µg

**Form:**

Available lyophilized

**Species:**

human

**Expression Source:**

*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**

Not tested

**Endotoxin Level:**

Please contact the lab for more information

**Peptide Sequence:**

KSTKLG~~LNQ~~NKKAAKLYKRERAQKRRNSSIPNGKLPME  
SPAQNFTIKEKDLE~~LDTE~~EPPNYDEEMSGGIEFLASV  
TKDTASDSPA~~GIDNPV~~FSPDEALDR~~SL~~LARLPPWLSPG  
ETVVPSQRARTQIPYSPGTFRR~~LMPF~~RLSSKSVDSFLQ  
ADGPEERPPAALPESTHM  
(665-834 aa encoded by BC101669)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.