

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human HIST3H2BB protein



Basic Information

Catalog Number:

Ag26901

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MPDPSKSA P A P K K G S K K A V T K A Q K K D G K K R K R G R K E S
Y S I V Y Y K V L K Q V H P D T G I S S K A M G I M N S F V N D I F E R I A
S E A S R L A H Y N K R S T I T S R E V Q T A V R L L L P G E L A K H A V S
E G T K A V T K Y T S S K

(1-126 aa encoded by BC100855)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.