

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human URAT1 protein

## Basic Information

**Catalog Number:**  
Ag27805

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PET28a, with N-terminal 6<sup>th</sup>His.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
SELLDLVGGGLGRFQVLQTMALMVSIMWLCTQSMLEN  
FSAAVPSHRCWAPLLDNSTAQASILGSLPEALLAISI  
PPGPNQRPHQCRRFRQPQWQLDPNATATSWSEAD  
TEPCVDGWVYDRSIFTSTIVAKWNLVCDSHALKPMA  
QSIYLAGILVGAAACGPASDRFGRRLLVTSYLQMAV  
MGTAARAFAPFPVYCLFRFLAFAVAGVMMNTGTL  
MEWTAARARPLVMTNLNLSLGFSGHGLTAAVAYGVRD  
WTLQLVVSVPFFLCFLYSWWLAESARWLLTTGRLD  
WGLQELWRVAIINGKGAVDTLTPEVLLSAMREELS  
MGQPPASLGTLLRMPGLRFRTC  
(4-351 aa encoded by BC053348)

## Reconstitution and Storage

**Reconstitution:**  
Reconstitute at 0.25 µg/ µl in 200 µl sterile water for short-term storage.  
After reconstitution with sterile water, if glycerol has no effect on subsequent experiments, it is recommended to add an equal volume of glycerol for long-term storage (see Stability and Storage for more details).  
If a different concentration is needed for your purposes please adjust the reconstitution volume as required (please note: the ion concentration of the final solution will vary according to the volume used).  
Note: Centrifuge vial before opening. When reconstituting, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution.

**Shipping:**  
The product is shipped at ambient temperature. Upon receipt, store it immediately at the recommended temperature (see below).

## Purity

85%, by SDS-PAGE with Coomassie Brilliant Blue staining.

## Formulation

The purified protein was Lyophilized from sterile PBS (58mM Na<sub>2</sub>HPO<sub>4</sub>, 17mM NaH<sub>2</sub>PO<sub>4</sub>, 68mM NaCl, pH7.4). 5 % trehalose and 5 % mannitol are added as protectant before lyophilization. The elution buffer contain 300mM imidazole.

## Stability and Storage


Store for up to 12 months at -20°C to -80°C as lyophilized powder.

## Storage of Reconstituted Protein

**Short Term Storage:**  
Store at 2-8°C for (1-2 weeks).

**Long Term Storage:**  
Aliquot and store at -20°C to -80°C for up to 3 months, buffer containing 50% glycerol is recommended for reconstitution. Avoid repeat freeze-thaw cycles.

## Selected Validation Data



For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.