

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human MATR3 protein



Basic Information

Catalog Number:

Ag31502

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, pET28a-SUMO, with N-terminal his.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

GRVIHLSNLPHSGYSDSAVLKLAEPYGIKKNYILMRMKS
QAFIEMETREDAMAMVDHCLKKALWFQGRGVKVDL
SEKYKKLVLRIPNRGIDLLKKDKSRKRSYSPDGKESPSDK
KSKTDGSQKTESSTEGKEQEEKSGEDGEKDTKDDQTE
QEPNMLLESEDELLVDEEEAALLESGLSSVGDETDLA
NLGDAVDGKKEPSDKAVIKDGSASAAKKKKKKVD
KIEELDQENEALENGIKNEENTEPGAESSENADDPN
KDTSENADGQSDENKDDYTIPDEYRIGPYQPNVPVGI
DYVIPKTGFYCKLCSLFYTNEEVAKNTHCSSLPHYQKL
KKFLNKLAEERRQKKET

(495-847 aa encoded by BC015031)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.