

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant mouse P62,Sqstm1 protein



Basic Information

Catalog Number:

Ag35064

Size:

50 µg

Form:

Available lyophilized

Species:

mouse

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MASFTVKA YLLGKEEATREIRRF SFCFSPEPEAEQAQA
AGPGPCERLLSRVAVLFP TLRPGGFQAHYRDEDEGDLV
AFSSDEELTMAMSVVKDDIFRIYIKEKKECRREHRPPC
AQEAPRNMVHPNVICDGCNGPVVGTRYKCSVCPDY
DLCSVCEGKGLHREHSLIFPNPFGHLSDFSFSRWRWR
KLLKHGHFGWPGWEMGPPGNWSPRPPRAGDGRPCPT
AESASAPPEDPNVNF LKNVGESVAAAALSP LGIEVDID
VEHGKRSRLTPTTPESSTGTEDKSNTQPSSCSSEVS
KPDGAGEGPAQSLTEQMKKIALESVGGQPEEQMESGN
CSGGDDDWTHLSKEVDPSTGELQSLQMPSEGPSS
LDPSQEGPTGLKEAALYPHLPPEADPRLIESLSQMLSM
GFSDEGGWLRLLQTKNYDIGAALDTIQYSKHPPPL
(1-442 aa encoded by NM_011018)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.