

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human NARFL protein



## Basic Information

Catalog Number:

Ag4217

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

QLNPTDTARKLTSFFKKIGVHFVFDTAFSRHFSLLESQ  
REFVRRFRGQADCRQALPLLASACPGWICYAEKTHG  
SFILPHISTARSPQQVMGSLVKDFFAQQQLTPDKIY  
HVTVMPCYDKLEASRPDFFNQEHQTRDVCVLTG  
EVFRLLEEEGVSLPDLEPAPLDSLCSGASAEPTSHRG  
GGSGGYLEHVFRAARELFGIHVAEVTYKPLRNKDFQ  
EVTLEKEGQVLLHFAMAYGFRNIQNLVQRLKRGRCPY  
HYVEVMACPSGCLNGGGQLQAPDRPSRELLQHVERL  
YGMVRAEAPEDAPGVQELYTHWLQGT DSECAGRLLH  
TQYHAVEKASTGLGIRW

(129-476 aa encoded by BC030248)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.