

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human USP46 protein



## Basic Information

**Catalog Number:**  
Ag4330

**Size:**  
50 µg

**Form:**  
Available lyophilized

**Species:**  
human

**Expression Source:**  
*e coli*-derived, PGEX-4T, with N-terminal GST.

**Biological Activity:**  
Not tested

**Endotoxin Level:**  
Please contact the lab for more information

**Peptide Sequence:**  
MTVRNIASICNMGTNASALEKDIGPEQFPINEHYFGL  
VNFNGNTCYCNSVLQALYFCRPFRENVLAYKAQQKKK  
ENLLTCLADLFHSIATQKKKVGVIIPPKFISRLRKENDL  
FDNYMQQDAHEFLNYLLNTIADILQEEKKQEKQNGKL  
KNGNMNEPAENNKPELTWVHEIFQGTLTNETRCLNC  
ETVSSKDEDFLDLSVDVEQNTSITHCLRDFSNTETLCS  
EQKYCYETCCSKQEAQKRMVKKLPMILALHLKRFKY  
MEQLHRYTKLSYRVVFPLELRLFNTSSDAVNLD RMYDL  
VAVVVHCGSGPNRGHYITIVKSHGFWLLFDDDI VEKI  
DAQAIEEFYGLTSDISKNSESGYILFYQSRE  
(1-366 aa encoded by BC037574)

## Selected Validation Data

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.