

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human ARHGEF7 protein



## Basic Information

Catalog Number:

Ag5344

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PET28a, with N-terminal 6\*His.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

KELELQILTEAIRNWEQDDIKTLGNVTYMSQVLIQCAG  
SEEKNERYLLFPNVLLMSASPRMSGFIYQGKLPPTG  
MTITKLESENHRNAFEISGSMIERILVSCNNQDLQE  
WVEHLQKQTKVTSVGNPTIKPHSVPSHTLPSHPVTPS  
SKHADSKPAPLTPAYHTLPHPSHHGTPHTTINWGPLE  
PPKTPKPWLSLCLRPAPPLRPSAALCYKEDLSKSPKTM  
KKLLPKRKPERKPSDEEFASRKSTAALEEDAQILKVEA  
YCTSAKTRQTLNSTWQGTDLMHNHVLAADDQPSLD  
SLGRRSSLSRLEPSDLSDESDYDSIWTAHSYRMGSTSR  
KSCCSYISHQN

(455-803aa encoded by BC060776)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.