

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human UGT1A10 protein



Basic Information

Catalog Number:
Ag6052

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PKG, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MARAGWTSPLVCVCLLLTCGFAEAGKLLVVPMDGS
HWFTMQSVVEKILRGHEVVVMPEVSWQLERSLNC
TVKTYSTSYTLEDQNRREFMVFAHAQWKAQAQSIFSL
LMSSSSGFLDLFFSHCRSLFNDRKLVEYLKESSFDAVFL
DPFDTCGLIVAKYFSLPSVVFTRGIFCHHLEEGAQCP
APLSYVPNDLLGFSDAMTFKERVVNHIVHLEDHLFCQ
YLFRNALEIASEILQTPVTA YDLYSHTSIWLLRTDFVLD
YPKPVMNMFIFGGINCHQGKPLPMEFEAYINASGEH
GIVVFSLGSVMSEIPEKKAMAIALDALGKIPQTVLWRYT
GTRPSNLANNTILVK
(1-350 aa encoded by BC053576)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.