

FOR IN VITRO RESEARCH USE ONLY.  
NOT FOR USE IN HUMANS OR ANIMALS.

# Recombinant human APOBEC3B protein



## Basic Information

Catalog Number:

Ag6115

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

*e coli*-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

MNPQIRNPMERMYRDTFYDNFENEPILYGRSYTWLCY  
EVKIKRGRSNLLWDTGVFRGQVYFEPQYHAEMCFLS  
WFCGNQLPAYKCFQITWVSWTPCPCVAKLAEFLS  
EHPNVTLTISAARLYYYWERDYRRALCRLSQAGARVKI  
MDYEEFAYCWENFVYNEGQQFMPWYKFDENYAFH  
RTLKEILRLRIFSVAFTAAMRSCASWTWFLCSWTRPR  
STGSLGSSPGAPASPGAVPGKCVRSFRRTHT  
(1-251 aa encoded by BC053859)

## Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.