

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human MTA3 protein



Basic Information

Catalog Number:
Ag6400

Size:
50 µg

Form:
Available lyophilized

Species:
human

Expression Source:
e coli-derived, PKG, with N-terminal GST.

Biological Activity:
Not tested

Endotoxin Level:
Please contact the lab for more information

Peptide Sequence:
MAANMYRVGDYVYFENSSNPYLIRRIEELNKTASGN
VEAKVVCFYRRRDISNTLIMLADKHAKIEEESETTVE
ADLTDKQKHQKHLKRELFSLRQYESLPATHIRGKCSVAL
LNETESVLSYLDKEDTFFYSLVYDPSLKTLLADKGEIRV
GPRYQADIPEMLLEGESDEREQSKLEVKVVDPNSPLT
DRQIDQFLVVARAVGTFFARALDCSSSVRQPSLHMSA
AAASRDITLFHAMDTLYRHSYDLSSAISVLVPLGGPVL
CRDEMEEWSASEASLFEAELEKYGKDFNDIRQDFLPW
KSLTSIIIEYYMWKTTDRYVQQKRLKAAEAESKLQVY
IPTYS
(1-343 aa encoded by BC053631)

Selected Validation Data

For technical support and original validation data for this product please contact:
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.