

For Research Use Only

CoraLite® Plus 488-conjugated TCOF1 Polyclonal antibody



Catalog Number: CL488-11003

Featured Product

Basic Information

Catalog Number:

CL488-11003

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1447

GenBank Accession Number:

BC014559

GeneID (NCBI):

6949

UNIPROT ID:

Q13428

Full Name:

Treacher Collins-Franceschetti
syndrome 1

Calculated MW:

97 kDa, 144 kDa, 152 kDa

Observed MW:

220 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : MCF-7 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite® Plus 488 TCOF1 antibody (CL488-11003) at dilution of 1:200.