

For Research Use Only

# CoraLite® Plus 488-conjugated PHLDA1 Polyclonal antibody



Catalog Number: CL488-18263

Featured Product

## Basic Information

Catalog Number:

CL488-18263

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG13125

GenBank Accession Number:

BC018929

GeneID (NCBI):

22822

UNIPROT ID:

Q8WV24

Full Name:

pleckstrin homology-like domain,

family A, member 1

Calculated MW:

45 kDa

Observed MW:

40-45 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : A375 cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926


E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

---

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed A375 cells using Coralite® Plus 488 PHLDA1 antibody (CL488-18263) at dilution of 1:200.