

For Research Use Only

CoraLite® Plus 488-conjugated CTTNBP2NL Polyclonal antibody



Catalog Number: **CL488-25523**

Basic Information

Catalog Number:

CL488-25523

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG22137

GenBank Accession Number:

BC016029

GeneID (NCBI):

55917

UNIPROT ID:

Q9P2B4

Full Name:

CTTNBP2 N-terminal like

Calculated MW:

639 aa, 70 kDa

Observed MW:

70-75 kDa

Purification Method:

Antigen affinity purification

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC (Intra)

Species Specificity:

human, mouse

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1x10⁶ MCF-7 cells were intracellularly stained with 0.8 ug Coralite® Plus 488 Cttbp2nl Polyclonal Antibody (CL488-25523)(red), or 0.8 ug CoraLite® Plus 488-conjugated Rabbit IgG control Rabbit PolyAb (CL488-30000) (blue). Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).