

For Research Use Only

CoraLite® Plus 488-conjugated IL-36 Beta/IL-1F8 Monoclonal antibody



Catalog Number: CL488-60290

Basic Information

Catalog Number: CL488-60290	GenBank Accession Number: BC101833	Purification Method: Protein G purification
Size: 1000 µg/ml	GeneID (NCBI): 27177	CloneNo.: 10D12B1
Source: Mouse	UNIPROT ID: Q9NZH7	Recommended Dilutions: IF-P 1:50-1:500
Isotype: IgG1	Full Name: interleukin 1 family, member 8 (eta)	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG12675	Calculated MW: 157 aa, 18 kDa	

Applications

Tested Applications:
IF-P
Species Specificity:
human

Positive Controls:
IF-P : human tonsillitis tissue,

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 IL-36 Beta/IL-1F8 antibody (CL488-60290, Clone: 10D12B1) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite® Plus 488 IL-36 Beta/IL-1F8 antibody (CL488-60290, Clone: 10D12B1) at dilution of 1:200.