

For Research Use Only

# CoraLite® Plus 488 Anti-Mouse CD117 (2B8)



Catalog Number: CL488-65054

## Basic Information

<b>Catalog Number:</b> CL488-65054	<b>GenBank Accession Number:</b> BC052457	<b>Purification Method:</b> Affinity purification
<b>Size:</b> 100ug, 500 µg/ml	<b>GeneID (NCBI):</b> 16590	<b>CloneNo.:</b> 2B8
<b>Source:</b> Rat	<b>UNIPROT ID:</b> P05532	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Isotype:</b> IgG2b, kappa	<b>Full Name:</b> kit oncogene	

## Applications

**Tested Applications:**  
FC

**Species Specificity:**  
Mouse

## Background Information

CD117, also known as c-Kit and SCFR, is a transmembrane protein with tyrosine kinase activity encoded by the oncogene c-kit (PMID: 2448137). It is a member of the type III receptor tyrosine kinase family, which also includes CSF-1R, PDGFR $\beta$ , PDGFR $\alpha$ , and FLT3 (PMID: 29518044). CD117 is expressed on hematopoietic stem cells and progenitor cells, mast cells, and is also found in a wide range of non-haemopoietic cell types (including melanocytes, germ cells, astrocytes, renal tubules, breast glandular epithelial cells, sweat glands, and interstitial cells of Cajal) (PMID: 10582338; 23073628). CD117 plays an important role in early haemopoiesis. It is also involved in pigmentation, fertility, gut movement, and some aspects of the nervous system (PMID: 23073628).

## Storage

**Storage:**  
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

**Storage Buffer:**  
PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1X10<sup>6</sup> mouse bone marrow cells were surface stained with APC Anti-Mouse CD45R (B220) (APC-65139, Clone: RA3-6B2) and CoraLite® Plus 488-conjugated Rat IgG2b isotype control (left) or 1 ug CoraLite® Plus 488 Anti-Mouse CD117 (CL488-65054, Clone: 2B8) (right). Cells were not fixed.