

For Research Use Only

CoraLite® Plus 488 Anti-Human CD1c (L161)



Catalog Number: CL488-65176

Basic Information

Catalog Number:

CL488-65176

Size:

100 tests, 5 µl/test

Source:

Mouse

Isotype:

IgG1, kappa

GenBank Accession Number:

BC126465

GeneID (NCBI):

911

UNIPROT ID:

P29017

Full Name:

CD1c molecule

Purification Method:

Affinity purification

CloneNo.:

L161

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

FC

Species Specificity:

Human

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide and 0.5% BSA.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 human PBMCs were surface stained with CoraLite® Plus 750 Anti-Human CD19 and 5 ul CoraLite® Plus 488 Anti-Human CD1c (L161) (CL488-65176) or CoraLite® Plus 488 Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.



1×10^6 human PBMCs were surface stained with 5 ul CoraLite® Plus 488 Anti-Human CD1c (L161) (CL488-65176) or CoraLite® Plus 488 Mouse IgG1 Isotype Control. Cells were not fixed. Lymphocytes were gated.