

For Research Use Only

# CoraLite® Plus 488 Anti-Mouse GL7 (GL-7)



Catalog Number: **CL488-65261**

## Basic Information

Catalog Number:

CL488-65261

Size:

100ug, 500 µg/ml

Source:

Rat

Isotype:

IgM, kappa

GenBank Accession Number:

GeneID (NCBI):

104230

Full Name:

Lymphocyte antigen 77

Purification Method:

Affinity purification

CloneNo.:

GL-7

Excitation/Emission maxima  
wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

FC

Species Specificity:

Mouse

## Background Information

### Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



1x10<sup>6</sup> unstimulated or anti-CD3/CD28 stimulated (3 days) Balb/c mouse splenocytes were surface co-stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) and 0.5 ug CoraLite® Plus 488 Anti-Mouse GL7 (CL488-65261, Clone:GL-7). Cells were not fixed.