

For Research Use Only

CoraLite® Plus 488-conjugated TrkC Monoclonal antibody



Catalog Number: **CL488-66380**

Basic Information

Catalog Number:

CL488-66380

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG2635

GenBank Accession Number:

BC013693

GeneID (NCBI):

4916

UNIPROT ID:

Q16288

Full Name:

**neurotrophic tyrosine kinase,
receptor, type 3**

Calculated MW:

839 aa, 94 kDa

Observed MW:

110 kDa, 160 kDa

Purification Method:

Protein G purification

CloneNo.:

2A7F11

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

Species Specificity:

human, mouse, rat, pig

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

