

For Research Use Only

CoraLite® Plus 488-conjugated DDX54 Monoclonal antibody



Catalog Number:CL488-66664

Basic Information

Catalog Number:

CL488-66664

Size:

1000 ug/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG25289

GenBank Accession Number:

BC156669

GeneID (NCBI):

79039

UNIPROT ID:

Q8TDD1

Full Name:

DEAD (Asp-Glu-Ala-Asp) box
polypeptide 54

Calculated MW:

98 kDa

Observed MW:

98 kDa

Purification Method:

Protein A purification

CloneNo.:

1H9F4

Recommended Dilutions:

IF/ICC 1:200-1:800

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC, FC (Intra)

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : PC-3 cells,

Background Information

DDX54, also named as ATP-dependent RNA helicase DDX54, is a 881 amino acid protein, which contains 1 helicase ATP-binding domain and belongs to the DEAD box helicase family. DDX54/DBP10 subfamily. DDX54 localizes in nucleus and Interacts in a hormone-dependent manner with nuclear receptors. DDX54 has RNA-dependent ATPase activity and represses the transcriptional activity of nuclear receptors.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed PC-3 cells using Coralite® Plus 488 DDX54 antibody (CL488-66664, Clone: 1H9F4) at dilution of 1:400.



1x10⁶ HepG2 cells were intracellularly stained with 0.4 ug Coralite® Plus 488-conjugated DDX54 Monoclonal antibody (CL488-66664, Clone:1H9F4) (red), or 0.4 ug Isotype Control (blue). Cells were fixed and permeabilized with Transcription Factor Staining Buffer Kit (PF00011).