

For Research Use Only

# CoraLite® Plus 488-conjugated NRBP2 Monoclonal antibody



Catalog Number: CL488-66687

## Basic Information

<b>Catalog Number:</b> CL488-66687	<b>GenBank Accession Number:</b> BC113873	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 340371	<b>CloneNo.:</b> 4F10H7
<b>Source:</b> Mouse	<b>UNIPROT ID:</b> Q9NSYO	<b>Recommended Dilutions:</b> IF/ICC 1:50-1:500
<b>Isotype:</b> IgG2a	<b>Full Name:</b> nuclear receptor binding protein 2	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Immunogen Catalog Number:</b> AG17011	<b>Calculated MW:</b> 501 aa, 58 kDa	
	<b>Observed MW:</b> 58 kDa, 33 kDa	

## Applications

<b>Tested Applications:</b> IF/ICC	<b>Positive Controls:</b> IF/ICC : HeLa cells, HEK-293 cells
<b>Species Specificity:</b> Human, Mouse, Rat, Pig	

## Background Information

Nuclear receptor binding protein 2 (NRBP2), also named as TRG16, belongs to the Ser/Thr protein kinase family. NRBP2 can regulate apoptosis of neural progenitor cells during their differentiation. Expression of NRBP2 was up-regulated following differentiation in rodent neural stem/progenitor cells. It has 4 isoforms with MW of 58, 33, 36 and 14 kDa. Endogenous NRBP2 in HEK293 cells has an apparent molecular mass of 55-60 kD by SDS-PAGE.

## Storage

**Storage:**  
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using CoraLite® Plus 488-conjugated NRBP2 antibody (CL488-66687, Clone: 4F10H7) at dilution of 1:100.



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CL488-66687 (NRBP2 antibody) at dilution of 1:100.