

For Research Use Only

# CoraLite® Plus 488-conjugated N-PAC Monoclonal antibody



Catalog Number: CL488-67340

## Basic Information

Catalog Number:

CL488-67340

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2b

Immunogen Catalog Number:

AG7133

GenBank Accession Number:

BC064940

GeneID (NCBI):

84656

UNIPROT ID:

Q49A26

Full Name:

cytokine-like nuclear factor n-pac

Calculated MW:

61 kDa

Observed MW:

60 kDa

Purification Method:

Protein A purification

CloneNo.:

1A5E1

Recommended Dilutions:

IF 1:300-1:1200

Excitation/Emission maxima wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

Human, Mouse, Rat

Positive Controls:

IF : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 N-PAC antibody (CL488-67340, Clone: 1A5E1) at dilution of 1:600.