

For Research Use Only

CoraLite® Plus 488-conjugated Phospho-TAU (Thr231) Recombinant antibody

Catalog Number: CL488-80101-7



Basic Information

Catalog Number: CL488-80101-7	GenBank Accession Number: BC000558	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 4137	CloneNo.: 240632B1
Source: Rabbit	UNIPROT ID: P10636	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Isotype: IgG	Full Name: microtubule-associated protein tau	
	Calculated MW: 37-46, 79-81 kDa	
	Observed MW: 50-80 kDa	

Applications

Tested Applications:
FC (Intra)

Species Specificity:
human

Background Information

Tau (tubulin-associated unit) is a microtubule-associated protein (also known as MAPT), expressed mainly in neurons of the central nervous system. Its primary function is to modulate microtubule dynamics for maintaining axonal cytoskeleton. The Tau protein has six isoforms produced from a single gene through alternative RNA splicing. Isoforms differ in number of inserts at the N-terminal half and the number of repeats at the C-terminal half (3 repeat-3R; 4 repeat-4R). Tau is hyperphosphorylated during aging and in age-related neurodegenerative diseases such as Alzheimer's disease (AD) and fronto-temporal dementia. Hyperphosphorylation of Tau leads to the formation of neurofibrillary tangles (NFT) in the neurons and glia cells, which is one of the hallmarks of AD.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ SH-SY5Y cells untreated (dashed lines) or treated with λ phosphatase which intracellularly stained with 0.25 μ g CoraLite[®] Plus 488 Phospho-TAU (Thr231) Recombinant Antibody (CL488-80101-7, Clone:240632B1)(red), or 0.25 μ g CoraLite[®] Plus 488 Rabbit IgG Isotype Control RecAb (CL488-98136, Clone: 240953C9) (blue). Cells were fixed with 4% PFA and permeabilized with 90% MeOH.