

For Research Use Only

CoraLite® Plus 488-conjugated HNRPDL Recombinant antibody

Catalog Number:CL488-83129-5



Basic Information

Catalog Number:

CL488-83129-5

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG1046

GenBank Accession Number:

BC007392

GeneID (NCBI):

9987

UNIPROT ID:

O14979

Full Name:

heterogeneous nuclear
ribonucleoprotein D-like

Calculated MW:

46 aa, 4 kDa

Observed MW:

38 kDa

Purification Method:

Protein A purification

CloneNo.:

230415F9

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima
wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse

Positive Controls:

IF/ICC : U-251 cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed U-251 cells using CoraLite® Plus 488 HNRPDL antibody (CL488-83129-5, Clone: 230415F9) at dilution of 1:200.