

For Research Use Only

# CoraLite® Plus 488-conjugated USP8 Recombinant antibody

Catalog Number:CL488-83172-2



## Basic Information

Catalog Number:

CL488-83172-2

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG27143

GenBank Accession Number:

BC110590

GeneID (NCBI):

9101

UNIPROT ID:

P40818

Full Name:

ubiquitin specific peptidase 8

Calculated MW:

1118 aa, 128 kDa

Purification Method:

Protein A purification

CloneNo.:

240016D3

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

## Applications

Tested Applications:

IF/ICC

Species Specificity:

human

Positive Controls:

IF/ICC : HeLa cells,

## Background Information

### Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HeLa cells using CoraLite® Plus 488 USP8 antibody (CL488-83172-2, Clone: 240016D3) at dilution of 1:200, CL594-Phalloidin (red).