

For Research Use Only

CoraLite® Plus 488-conjugated THOC5 Recombinant antibody

Catalog Number:CL488-83435



Basic Information

Catalog Number: CL488-83435	GenBank Accession Number: BC003615	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 8563	CloneNo.: 240013F7
Source: Rabbit	UNIPROT ID: Q13769	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG	Full Name: THO complex 5	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG6664	Calculated MW: 79 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : U2OS cells,
Species Specificity: human	

Background Information

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Acetone) fixed U2OS cells using CoraLite® Plus 488 THOC5 antibody (CL488-83435, Clone: 240013F7) at dilution of 1:200.