

For Research Use Only

CoraLite® Plus 488-conjugated KMT2C Recombinant antibody

Catalog Number:CL488-83595-4



Basic Information

Catalog Number: CL488-83595-4	GenBank Accession Number: NM_170606	Purification Method: Protein A purification
Size: 1000 µg/ml	GeneID (NCBI): 58508	CloneNo.: 240542A2
Source: Rabbit	UNIPROT ID: Q8NEZ4	Recommended Dilutions: IF/ICC 1:50-1:500
Isotype: IgG	Full Name: myeloid/lymphoid or mixed-lineage leukemia 3	Excitation/Emission maxima wavelengths: 493 nm / 522 nm
Immunogen Catalog Number: AG28875	Calculated MW: 541 kDa	

Applications

Tested Applications: IF/ICC	Positive Controls: IF/ICC : MCF-7 cells,
Species Specificity: human	

Background Information

Lysine (K)-specific methyltransferase 2C (KMT2C, also known as MLL3) belongs to the mixed-lineage leukemia (MLL) family of histone methyltransferases which methylate the histone 3 tail at lysine 4 (H3K4) as part of the complex proteins associated with Set 1 (COMPASS) complex. Although originally identified as oncogenic fusions in leukemia, genome-wide mutation studies have revealed frequent, presumably loss-of-function, mutations in various members of the MLL family, including MLL2/KMT2B, MLL3/KMT2C, and MLL4/KMT2D in a variety of malignancies, particularly solid tumors. Mechanistic studies of KMT2C in normal cells have focused primarily on its role in enhancer regulation by deposition of H3K4me1 marks.

Storage

Storage:
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.
Storage Buffer:
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed MCF-7 cells using CoraLite® Plus 488 KMT2C antibody (CL488-83595-4, Clone: 240542A2) at dilution of 1:200.