

For Research Use Only

# CoraLite® Plus 488-conjugated DLAT Recombinant antibody



Catalog Number: CL488-83654-3

## Basic Information

<b>Catalog Number:</b> CL488-83654-3	<b>GenBank Accession Number:</b> BC039084	<b>Purification Method:</b> Protein A purification
<b>Size:</b> 1000 µg/ml	<b>GeneID (NCBI):</b> 1737	<b>CloneNo.:</b> 240713C12
<b>Source:</b> Rabbit	<b>UNIPROT ID:</b> P10515	<b>Recommended Dilutions:</b> IF/ICC 1:50-1:500
<b>Isotype:</b> IgG	<b>Full Name:</b> dihydrolipoamide S-acetyltransferase	<b>Excitation/Emission maxima wavelengths:</b> 493 nm / 522 nm
<b>Immunogen Catalog Number:</b> AG22831	<b>Calculated MW:</b> 647 aa, 70 kDa	
	<b>Observed MW:</b> 69-70 kDa	

## Applications

<b>Tested Applications:</b> IF/ICC	<b>Positive Controls:</b> IF/ICC : HepG2 cells,
<b>Species Specificity:</b> human, mouse, rat	

## Background Information

### Storage

**Storage:**  
Store at -20°C. Avoid exposure to light. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using Coralite® Plus 488 DLAT antibody (CL488-83654-3, Clone: 240713C12 ) at dilution of 1:200.