

For Research Use Only

CoraLite® Plus 488-conjugated HINT2 Recombinant antibody

Catalog Number:CL488-84031-5



Basic Information

Catalog Number:

CL488-84031-5

Size:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG12413

GenBank Accession Number:

BC047737

GeneID (NCBI):

84681

UNIPROT ID:

Q9BX68

Full Name:

histidine triad nucleotide binding

protein 2

Calculated MW:

163 aa, 17 kDa

Observed MW:

14 kDa

Purification Method:

Protein A purification

CloneNo.:

241150D10

Recommended Dilutions:

IF/ICC 1:50-1:500

Excitation/Emission maxima

wavelengths:

493 nm / 522 nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

human, mouse, rat

Positive Controls:

IF/ICC : HeLa cells,

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed HeLa cells using CoraLite® Plus 488 HINT2 antibody (CL488-84031-5, Clone: 241150D10) at dilution of 1:200, CL594-phalloidin (red).