

For Research Use Only

CoraLite® Plus 555 Anti-Mouse GL7 (GL-7)



Catalog Number: **CL555-65261**

Basic Information

Catalog Number:

CL555-65261

Size:

100ug, 0.5 mg/ml

Source:

Rat

Isotype:

IgM, kappa

GenBank Accession Number:

GeneID (NCBI):

104230

Full Name:

Lymphocyte antigen 77

Purification Method:

Affinity purification

CloneNo.:

GL-7

Excitation/Emission maxima
wavelengths:

554 nm / 570 nm

Applications

Tested Applications:

FC

Species Specificity:

Mouse

Background Information

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ unstimulated or anti-CD3/CD28 stimulated mouse splenocytes were surface co-stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) and 0.5 ug CoraLite® Plus 555 Anti-Mouse GL7 (CL555-65261, Clone:GL-7). Cells were not fixed.