

For Research Use Only

CoraLite®555-conjugated VPS53 Monoclonal antibody



Catalog Number:CL555-67610

Basic Information

Catalog Number:

CL555-67610

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG29875

GenBank Accession Number:

BC029560

GeneID (NCBI):

55275

UNIPROT ID:

Q5VIR6

Full Name:

vacuolar protein sorting 53 homolog
(S. cerevisiae)

Calculated MW:

832 aa, 94 kDa

Observed MW:

94 kDa

Purification Method:

Protein G purification

CloneNo.:

2B12C3

Recommended Dilutions:

IF/ICC 1:200-1:800

Excitation/Emission maxima
wavelengths:

557 nm / 570nm

Applications

Tested Applications:

IF/ICC

Species Specificity:

Human, Mouse, Rat

Positive Controls:

IF/ICC : HEK-293 cells,

Background Information

VPS53 is a component of the Golgi-associated retrograde protein (GARP) complex, also called VFT (VPS fifty-three) complex, composed of VPS51, VPS52, VPS53 and VPS54. GARP complex functions in traffic from endosomes to the trans-Golgi network (PMID: 15878329). GARP proteins interact with RAB proteins and SNARE proteins.

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (-20°C Ethanol) fixed HEK-293 cells using CoraLite@555 VPS53 antibody (CL555-67610, Clone: 2B12C3) at dilution of 1:400.