

For Research Use Only

CoraLite®594 Anti-Mouse CD44 (IM7)



Catalog Number: **CL594-65117**

Basic Information

Catalog Number: CL594-65117	GenBank Accession Number: BC051388	Purification Method: Affinity purification
Size: 100ug, 500 µg/ml	GeneID (NCBI): 12505	CloneNo.: IM7
Source: Rat	UNIPROT ID: P15379	Excitation/Emission maxima wavelengths: 588 nm / 604 nm
Isotype: IgG2b, kappa	Full Name: CD44 antigen	

Applications

Tested Applications:
FC

Species Specificity:
Mouse

Background Information

CD44 is a type I transmembrane glycoprotein expressed on embryonic stem cells and in various levels on other cell types including connective tissues and bone marrow. CD44 expression is also upregulated in subpopulations of cancer cells and is recognized as a molecular marker for cancer stem cells (PMID: 29747682). It is a cell-surface receptor that mediates cell-cell and cell-matrix interactions through its affinity for hyaluronic acid (HA) and possibly also through its affinity for other ligands (PMID: 10694938). Adhesion with HA plays an important role in cell migration, tumor growth and progression. CD44 is also involved in lymphocyte activation, recirculation and homing, and in hematopoiesis.

Storage

Storage:
Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:
PBS with 0.1% sodium azide and 0.5% BSA, pH 7.3.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface co-stained with CoraLite® Plus 647 Anti-Mouse CD3 (17A2) and 0.5 ug CoraLite®594 Anti-Mouse CD44 (CL594-65117, Clone: IM7) or 0.5 ug CoraLite®594 Rat IgG2b Isotype Control (LTF-2) (CL594-65211, Clone: LTF-2). Cells were not fixed.