

For Research Use Only

CoraLite®594-conjugated SLC9A9 Monoclonal antibody



Catalog Number:CL594-66577

Basic Information

Catalog Number:

CL594-66577

Size:

1000 µg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG4660

GenBank Accession Number:

BC035779

GeneID (NCBI):

285195

UNIPROT ID:

Q8IVB4

Full Name:

solute carrier family 9
(sodium/hydrogen exchanger),
member 9

Calculated MW:

645 aa, 73 kDa

Observed MW:

73 kDa

Purification Method:

Protein A purification

CloneNo.:

1B8D9

Recommended Dilutions:

WB 1:500-1:1000

Excitation/Emission maxima
wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

WB

Species Specificity:

Human, Mouse, Rat

Positive Controls:

WB : HeLa cells, HEK-293 cells, HepG2 cells

Background Information

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% Glycerol, 0.05% Proclin300, 0.5% BSA, pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with CL594-66577 (SLC9A9 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.