

For Research Use Only

CoraLite®594 Anti-Mouse CD86 Rabbit Recombinant Antibody

Catalog Number:CL594-98025



Basic Information

Catalog Number:

CL594-98025

Size:

100ug, 500 µg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

NM_019388

GeneID (NCBI):

12524

UNIPROT ID:

P42082-1

Full Name:

CD86 antigen

Calculated MW:

35 kDa

Purification Method:

Protein A purification

CloneNo.:

230476B8

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

FC

Species Specificity:

mouse

Background Information

CD86 (also known as B7-2) is a costimulatory molecule belonging to the immunoglobulin superfamily. Primarily expressed on antigen-presenting cells (APCs), including B cells, dendritic cells, and macrophages, CD86 is the ligand for two proteins at the cell surface of T cells, CD28 antigen and CTLA-4. Binding of CD86 with CD28 antigen is a costimulatory signal for activation of the T-cell. Binding of CD86 with CTLA-4 negatively regulates T-cell activation and diminishes the immune response. (PMID: 7513726; 1847722; 11029388; 27591335)

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide.

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1×10^6 untreated or LPS treated mouse splenocytes were surface stained with 0.25 ug CoraLite®594 Anti-Mouse CD86 Rabbit Recombinant Antibody (CL594-98025, Clone: 230476B8). Cells were not fixed.



1×10^6 LPS treated mouse splenocytes were surface stained with 0.25 ug CoraLite®594 Anti-Mouse CD86 Rabbit Recombinant Antibody (CL594-98025, Clone: 230476B8) (red) or 0.25 ug Isotype Control (blue). 1×10^6 untreated mouse splenocytes were surface stained with 0.25 ug CoraLite®594 Anti-Mouse CD86 Rabbit Recombinant Antibody (CL594-98025, Clone:230476B8) (black). Cells were not fixed.